



Australia's premium manufacturer of medical and laboratory products

Product Specification Sheet

Series Peltier Cooled Incubators

Model ip-80

80 litre peltier cooled incubator



*Optional double glazed door shown





Australia's premium manufacturer of medical and laboratory products

Product Specification Sheet

Series Peltier Cooled Incubators

Model ip-80

80 litre peltier cooled incubator

Chamber capacity 80 Litres
Chamber width 435 mm
Chamber depth 375 mm
Chamber height 502 mm

Chamber material Stainless Steel 304L

External width 541 mm
External depth 575 mm
External height 760 mm

External material Durable epoxy powder-coated steel

Recommended clearance 75mm on sides and rear

Temperature range (15C below ambient) to 50C

Temperature sensor Pt100
Control accuracy 0.1 C
Temporal stability 0.2 C

Available shelf positions 4
Shelves supplied 2

Shelf dimensions 428 x 320 mm

Shelf spacing 100mm Shelf capacity 15 KG

Insulation type Mineral Wool - asbestos free

Insulation thickness 50mm

Maximum cooling capacity (watts) 200 Maximum heating capacity (watts) 200

Steridium Pty Ltd | Brisbane | Australia | www.steridium.com | sales@steridium.com | +61 7 55 931 352





Australia's premium manufacturer of medical and laboratory products

Product Specification Sheet

Series Model Peltier Cooled Incubators

ip-80

80 litre peltier cooled incubator

Electrical supply Single phase, 220 - 240V, 50 Hz

(alternatives available on request)

Circulating fan 1

Countdown timer 99 hr 59 min - or continuous

Controller type Multifunction intelligent digital PID-microprocessor
Readout Digital LCD readout of set and actual temperature

Temperature alarms High | low | deviation

Country of origin Australia

Options / Accessories

- Additional shelves
- 8 Step programmable controller
- 500 Step programmable controller
- Inner viewing door perspex
- Inner viewing door toughened glass
- Double glazed glass door
- Side access port 50 mm
- Internal power supply*
- Inner viewing light LED*
- RS232 or RS485 comms module
- Ethernet comms module
- ISPM 15 compliant crating for export
- Calibration certificate

Steridium Pty Ltd | Brisbane | Australia | www.steridium.com | sales@steridium.com | +61 7 55 931 352